



Cavex Oral Pre Rinse

Evaluation of effectiveness.

GENERAL:

Mouthrinses are used in dentistry for a number of reasons: to freshen breath, to control tooth decay, to reduce plaque formation, to prevent or reduce gingivitis, or to deliver a combination of these effects. Given their antiseptic properties, however, some agents also can serve an important function in infection control. The use of a mouthrinse by the patient before dental procedures is based on a similar principle of reducing the number of oral microorganisms. This reduction also reduces the number of microorganisms that may escape a patient's mouth during dental care through aerosols, spatter, or direct contact.

HYDROGEN PEROXIDE:

Systematic reviews to evaluate the efficacy of preprocedural mouthrinses have been performed and the reduction of the number of microorganisms released by in the aerosol are compared. Mouthrinses with chlorhexidine, essential oils, and cetylpyridinium chloride reduced the number of microorganism with more then 60%. However, due to the fact that corona type viruses are vulnerable to oxidation, preprocedural mouthrinse containing oxidative agents such as 1% hydrogen peroxide are recommended. A preprocedural mouthrinse would be even more useful in cases when rubber dam cannot be used.

The effectiveness of hydrogen peroxide against microorganisms has been proven for many years. As a surface disinfectant or in vapor form it is used in food, veterinary and healthcare. Hydrogen peroxide works by producing destructive radicals that can attack essential cell components (i.e. membrane lipids, DNA, and other essential cell components) As the degradation products are water and oxygen we can consider an aqueous peroxide product as environmental friendly.

Special interest for studies on the effective disinfection of institutions are performed, measuring the virucidal efficacy of peroxide against respiratory and enteric viruses. Hydrogen peroxide vapor shows high efficacy against the enveloped RNA viruses avian flu virus (AIV), swine influenza virus (SwIV) and the transmissible gastroenteritis corona virus (TGEV, a SARS-CoV surrogate) A solution of 0.5% hydrogen peroxide demonstrated already bactericidal and virucidal activity in 1 minute.

DENTAL GUIDELINES:

Although there are several dental guidelines published, showing differences in the details, the general information is comparable. Regarding patients pre rinse the guidelines are giving clear directions. In order to obtain virucidal effectiveness against the Covid-19 virus the only option is to pre rinse with 1 – 2 % hydrogen peroxide for 30 – 60 seconds.



ADDITIONS:

It has been reported that the ACE2 (angiotensin-converting enzyme 2) present in the oral mucosa is the main host cell receptor of Covid-19 and plays a crucial role in the entry of virus. Besides, it is proven that viruses can remain viable and infectious in aerosols for hours and on surfaces up to days. Therefore, the oral cavity is a potentially high risk and needs appropriate cleaning. Water based solutions of hydrogen peroxide create a thin foam in the mouth. By additions of a bio-adhesive the viscosity and adhesion of the foam can be optimized resulting in adherence especially to the glands in the mucosa, the source where virus containing saliva is exposed.

Addition of erythritol will result in not only virucidal activity but also anti-bacterial effectiveness. A very small quantity of glycerin will overcome a dry mouth feeling and together with a touch of mint flavor it will increase the acceptance by patients to perform the full rinsing session.

STORAGE:

Under normal conditions, hydrogen peroxide is sufficient stable when properly stored, preferably refrigerated in dark containers.

SAFETY:

Peroxides decompose in water and oxygen while releasing active radicals. As a result such solution possess anti microbial and cleaning effects. High concentrations of peroxide and long contact times can be destructive for soft tissues and metals. 1,5 % hydrogen peroxide for oral pre procedural treatment with short contact times (30-60 seconds) can be considered as being acceptable when handled by dental professionals. Always avoid contact with eyes and do not swallow.

CONCLUSION:

Cavex Oral Pre Rinse is developed to be used as a pre procedural rinse before starting the actual dental procedure. The 1,5 % hydrogen peroxide works as a bleaching agent, cleaning the oral cavity and reducing the total amount of microorganism. The addition of an bio adhesive in order to form the ideal foam condition and erythritol with anti bacterial properties will strengthen the cleaning effect. Small quantity of glycerin against dry mouth and slight mint flavor optimizes the patient product acceptance.

Date: 20 April 2020

Richard Woortman
Manager Technical Services

**REFERENCES:****REF 01 - Virus inactivation by hydrogen peroxide**

Vopr Virusol 1977 Nov-Dec;(6):731-3.

Mentel' R, Shirmakher R, Kevich A, Drežin RS, Shmidt I.

REF 02 - Virucidal efficacy of hydrogen peroxide vapour disinfection

Journal of Hospital Infection - Volume 80, Issue 2, February 2012, Pages 110-115

E.Tuladhar, P.Terpstra, M.Koopmans E.Duizer

REF 03 – Evaluating the virucidal efficacy of hydrogen peroxide vapour

Journal of Hospital Infection - Volume 86, 2014, Pages 255-259

SM Goya, Y Chander, S Yezli, JA Otter

REF 04 - COVID-19: Hype? Hazard? Dental Office Implications?

Clinicians Report Volume 13 Issue 5 May 2020

REF 05 - ADA Interim Guidance for Management of Emergency and Urgent Dental Care

American Dental Association April 2020

REF 06 – Infectiepreventive in de mondzorg tijdens uitbraak SARS-CoV-2

Ned Tijdschr Tandheelkd 2020; 127

CMC Volgenant, IF Persoon, RAG de Ruijter, JJ de Soet

REF 07 High expression of ACE2 receptor of 2019-nCoV on the epithelial cells of oral mucosa

International Journal of Oral Science (2020) 12:8

Hao Xu, Liang Zhong, Jiaxin Deng, Jiakuan Peng, Hongxia Dan, Xin Zeng, Taiwen Li and Qianming Chen

REF 08 Aerosol and surface stability of HCoV-19 (SARS-CoV-2) compared to SARS-CoV-1

The New England Journal of Medicine

Neeltje van Doremalen, Trenton Bushmaker, Dylan H Morris, Myndi G Holbrook, Amandine Gamble, Brandi N Williamson, Azaibi Tamin, Jennifer L Harcourt, Natalie J Thornburg, Susan I Gerber, James O Lloyd-Smith, Emmie de Wit, Vincent J Munster

REF 09 Transmission routes of 2019-nCoV and controls in dental practice

International Journal of Oral Science 2020 12:9

Xian Peng, Xin Xu, Yuqing Li, Lei Cheng, Xuedong Zhou and Biao Ren

REF 10 Aerosols and splatter in dentistry

A brief review of the literature and infection control Implications

JADA, Vol. 135, April 2004

Stephen K. Harrel, John Molinari

REF 11 Efficacy of preprocedural mouthrinses in the reduction of microorganisms in aerosol. A systematic review.

JADA 2019;150(12):1015-1026

Vanessa Costa Marui, Maria Luisa Silveira Souto, Emanuel Silva Rovai, Giuseppe Alexandre Romito, Leandro Ch

REF 12 Persistence of coronaviruses on inanimate surfaces and their inactivation with biocidal agents.

Journal of Hospital Infection - Volume 104, 2020, Pages 246-251

G Kampf, D Todt, S Pfaender, E Steinmann

REF 13 Development and characterisation of a moisture-activated bioadhesive drug delivery system for percutaneous local anaesthesia.

International Journal of Pharmaceutics 1998

AD Woolfson, DF McCafferty, GP Moss

REF 14 Erythritol is more effective than Xylitol and Sorbitol in managing oral health endpoints.

Journal of Dentistry Vol 2016

P de Cock, K akinen, E Honkala, M Saag, E Kennepohl, A Eapen