



■ **Cavex Holland BV**
 Fustweg 5
 2031 CJ Haarlem
 The Netherlands

Tel +31 23 530 7700
 dental@cavex.nl
 www.cavex.nl

Tulip alginate

the best choice in economical alginates!

Do you want to make an initial impression to prepare for the prosthetic treatment? Or the mould for the prosthesis itself? Or documentation or working models? With Tulip you can work any way you please.



Dust-free dosing and mixing. Rapid water absorption gives a smooth, homogeneous alginate paste without lumps, in just a short time. The right basis and the ideal consistency for a detailed impression. The Tulip paste is thixotropic. A consistency that flows easily, giving clear reproduction of all the contours of the jaw without the material running out of the tray.

The impression is stable, elastic and is easily removed from the mouth despite undercuts. The alginate impression gives high detail reproduction

(50µm) and is clearly legible due to the yellow colour of the material. Its high gypsum compatibility guarantees a nice, smooth surface to the gypsum. **Tulip has all the features of a high-quality alginate. Conforms with the ISO and ADA standards.**

Application specifications

	Tulip	Tulip Colorswitch
mixing time	30 sec	30 sec
working time	1.30 min	1.30 min
setting time in mouth	1.30 min	1.00 min
total setting time	3 min	2.30 min
mixing ratio (powder/water)	23g / 50ml	23g / 50ml

Tulip Colorswitch

impression material with a build-in timer



Tulip Colorswitch has all the advantages of Tulip alginate. Plus one unique feature: the changing colour. This helps you to time the process of mixing and taking the impression without being distracted from treatment or patient. Dentists are discovering another useful tool. The colour change of a little ball of Tulip Colorswitch in the patient's hand provides distraction and relaxation, resulting in a better impression.



At the beginning of the mixing process, the material is still violet. When the colour changes to dark pink the paste is sufficiently mixed.



The paste is ready for loading the tray. Insert the tray with pink paste in the mouth.



Wait until the paste turns off-white. This is the time to remove the tray from the mouth.



The off-white impression turns pink again after 15 - 30 minutes, so the dental technician receives a clearly legible impression and can ensure it is poured out properly.