

Cavex Bite&White StainLess**DESCRIPTION**

Teeth Whitening Activator & Stain Removal Paste. Removes surface stains and plaque from teeth. Enhances whitening treatments by optimizing the pH in the mouth.

- Mildly abrasive to enamel / dentin
- High cleaning and polishing values.
- Increases mouth pH as pre-treatment for teeth whitening.
- Not intended for regular use.

PACKAGING

Tube 35 ml – 50 g

COMPOSITION

Aqua, Sorbitol, Hydrated silica, Glycerin, Erythritol, Lithium magnesium sodium silicate, Poloxamer, Sodium hexametaphosphate, Cellulose gum, Xanthan gum, Phenoxyethanol, Sodium saccharin, Mint aroma, Limonene, CI 42090

INDICATIONS FOR USE

Before in-office and/or at-home tooth whitening

CONTRA-INDICATIONS

Do not use in cases of allergies to mint aroma or any of the ingredients.

HANDLING PROCEDURE

1. Dispense Cavex Bite&White StainLess onto toothbrush
2. Brush the teeth thoroughly
3. Clean teeth at least 2-3 minutes
4. Spit out any excess material. **Do not rinse with water**

WARNINGS & FIRST AID INSTRUCTIONS

In case of contact with the eyes, flush thoroughly with plenty of water and seek medical attention.

Do not use the material after the expiry date.

STORAGE ADVICE

Store in a cool place, dry place (frost-free) and protected from direct sunlight.

Made in The Netherlands

Our technical advice, whether verbal, in writing or by way of trials, is given in good faith but without warranty, and this also applies where proprietary rights of third parties are involved. It does not release you from the obligation to test the products supplied by us as to their suitability for the intended processes and uses.

The application, use and processing of the products are beyond our control and, therefore, entirely your own responsibility. Should, in spite of this liability be established for any damage, it will be limited to the value of the goods delivered by us and used by you. We will, of course, provide consistent quality of our products within the scope of our General Conditions of Sale and Delivery.