



Good oral hygiene is essential to good health. To safeguard your health, Cavex offers effective solutions that promote health for everyone and products that seamlessly fit the needs of modern practices and contemporary patients around the world.

Your impression is our concern

A revolutionary new disinfecting and cleaning fluid for alginate, polyether and silicone impressions. Safe, fast, effective.

The perfect alginate and silicone disinfection fluid, effectively kills bacteria, fungi and viruses but does not affect the surface of the impression. Cavex ImpreSafe is a non-toxic disinfectant that is aldehyde free and needs no more than 3 minutes to be effective. Due to this short contact time the surface of the impression will remain 100% intact and the end result is not compromised. The Cavex ImpreSafe system consists of 1 liter of Cavex ImpreSafe, a disinfection container, shipping bags and a timer, to ease the use.



Cavex ImpreSafe; cleaner and disinfectant for alginate, polyether and silicone impressions.

- **Fast and efficient**

Cavex ImpreSafe only needs 3 minutes to kill bacteria, fungi and viruses. Due to this short contact time and its unique composition, contrary to other disinfection fluids **Cavex ImpreSafe has no negative effect on the impression surface.**

- **Concentrate**

Cavex ImpreSafe is a highly concentrated disinfection fluid, saving storage space. 30 ml of Cavex ImpreSafe is sufficient for 1 liter of disinfectant. **So 1 bottle of Cavex ImpreSafe gives 33 liters of ready-to-use disinfection fluid.** The bottle has an easy dosage system.

- **Safe**

Cavex ImpreSafe is 100% aldehyde free making this non-toxic fluid very safe to use. Cavex ImpreSafe has a friendly, neutral smell.

- **Complete system**

Cavex ImpreSafe is available as a complete system. This system comprises 1 liter of Cavex ImpreSafe disinfectant, a disinfection container, shipping bags, a timer and a protocol.

- **Reliable and secure**

Cavex ImpreSafe has been thoroughly tested by expert laboratories (official reports on www.cavex.nl) and has proven to be a reliable and highly effective disinfection system. **Complies with DGHM guidelines and EN 1040, EN 1275, EN 13727, EN 13624, EN 14561, EN 14562.**

Cavex ImpreSafe and accessories are available from your dental supplier. For further information contact Cavex, +31 23 530 77 00 or www.cavex.nl



Assortment:

- HA053** Cavex ImpreSafe Starterkit
- HA050** Cavex ImpreSafe Disinfection fluid
- HA020** Cavex ImpreSafe Container