

PREMIUM ORAL CARE


Cavex VacuFormer System Accessories

Cavex Mouth Protector laminates

Professional 3.8mm dual layer and 5mm triple layer laminates with a combined hard and soft layer.

Cavex Mouth Protectors (3.8mm, double layer)

colour	article no.	content
Red/White/Blue	VT020	12 pcs
Clear	VT036	12 pcs



Cavex Mouth Protectors X-Pro (5mm, triple layer)

colour	article no.	content
Green	VT031	12 pcs
Transparent	VT033	12 pcs

Cavex Mouth Protectors Multi Color (4mm)

colour	article no.	content
Pretty Pink	VT037	6 pcs
FreeStyle	VT038	6 pcs
Splash	VT039	6 pcs



Professional mouthguards with two layers of soft EVA. The material has excellent tensile strength.



Cavex Bleach laminates

Made of a soft/clear material that is easy to cut, resulting in a perfect finish. The soft material is very comfortable on the gum tissue and thus ideal for bleaching trays.

Cavex Bleach (1,0mm)

colour	article no.	content
Transparent	VT022	25 pcs

Cavex Bleach Heavy (1,5mm)

colour	article no.	content
Transparent	VT030	25 pcs



Cavex Bite&White Barrier/Spacer

An easy-to-use blue colored resin specially developed for creating reservoirs in bleaching trays.

Cavex Bite&White Barrier/Spacer

article no.	content
VT025	2 pcs



Cavex Splint laminates

Made of a transparent, rigid material ideally suited for splints, surgical trays, orthodontic retainers, etc.

Cavex Splint (0,6mm)

colour	article no.	content
Transparent	VT021	50 pcs

Cavex Splint + (1,0mm)

colour	article no.	content
Transparent	VT028	25 pcs



Cavex Splint X-heavy (1,5mm)

colour	article no.	content
Transparent	VT026	25 pcs

Cavex Bruxism laminates

Hard and crystal clear, 2mm thickness. Perfect for the treatment of bruxism and helps to prevent wear and tear of the elements.



Cavex Bruxism (2,0mm)

colour	article no.	content
Transparent	VT027	25 pcs

YOUR
IMPRESSION
IS OUR
SPECIALTY

CAVEX