

What are the dimensions of the different sheets?

	Splint	Bleach	Mouthguard	Bruxism Splint X Heavy
dimensions (inch)	5 x 5	5 x 5	5 x 5	5 x 5
dimensions (mm)	12.7 x 12.7	12.7 x 12.7	12.7 x 12.7	12.7 x 12.7
thickness (inch)	0.022\0.039\0.059	0.039\0.059	0.15\0.16\0.20	0.079
thickness (mm)	0.6\1\1.5	1\1.5	3.8\4\5	2
Composition	Acrylic	Soft EVA	EVA	Acrylic