

PREMIUM ALGINATE SYSTEM

CAVEX IMPRESAFE

-
Disinfection system for alginate,
silicone and polyether impressions.

-
Kills 99.99% of all bacteria, fungi
and viruses in just 3 minutes.



"I loved the whole system.
It was very easy to use and
implement into our office."

YOUR IMPRESSION IS OUR SPECIALTY

CAVEX

PREMIUM ALGINATE SYSTEM

CAVEX IMPRESAFE

The Cavex ImpreSafe System is based on a non-toxic and aldehyde-free disinfection concentrate with an added booster and surfactant. This unique liquid kills 99.99% of all bacteria, fungi and viruses on an alginate or silicone impression in just 3 minutes, without affecting the surface of the impression.

The Cavex ImpreSafe starter kit consists of a bottle of disinfection concentrate, a container with a removable tray, a timer and 50 zip-lock bags for the perfect storage of disinfected impressions. The container and liquid are also available separately.

Thanks to its unique composition, Cavex ImpreSafe - unlike other disinfecting liquids - is effective in a very short contact time, so it has no negative effect on the surface of the impression.

Cavex ImpreSafe comes as a concentrate in a 1 liter bottle with an integrated dosing system. Because 30 ml of concentrate is sufficient for 1 liter of disinfectant liquid with which the container is filled, one bottle is sufficient for 33 liters of disinfectant liquid. Cavex ImpreSafe is 100% aldehyde-free, so non-toxic and extremely safe to use. In addition, the liquid has a pleasant, neutral scent.



Cavex ImpreSafe has been extensively tested by several leading laboratories and has proven to be a reliable and highly effective disinfectant liquid. The product complies with the DGHM directives and the standards EN 1040, EN 1275, EN 1276, EN 1650, EN 13624, EN 13727, EN 14476, EN 14561 and EN 14562 (see cavex.nl/impresafe for the test reports).

Characteristics

- A complete, fast and effective system
- A safe and reliable disinfection solution
- No negative effect on the print



Clinical rating of 96% by The Dental Advisor

We value the opinion of the professional. That is why we let The Dental Advisor test our Cavex ImpreSafe. After a clinical evaluation by a team of professionals, this product received a clinical rating of no less than 96%, which corresponds to the maximum score of 5 pluses. Professionals praise its simplicity, effectiveness and pleasant fragrance.

YOUR IMPRESSION IS OUR SPECIALTY

Fustweg 5 | 2031 CJ Haarlem | The Netherlands | cavex.nl

CAVEX

ZYH003EN.02